

Caregiver® Thermometer Mounting and Security Options





WALL MOUNT MODEL WMT-TF

The Caregiver® Wall Mount provides a versatile solution for storage of thermometers in the PRO-TF Series. It may be secured to a wall, mounted on commonly used "rail" systems or attached to an IV pole basket with optional accessories. The system consists of custom molded wall mount, attachment pin for use with optional tethering devices, mounting screws and drywall anchors.



BASKET MOUNT ADAPTER MODEL BSM-TF

The Basket Mount Adapter is used with the Wall Mount WMT-TF. It permits attachment of the Wall Mount to baskets typically used on IV poles.



SECURE TETHER MODEL SMT-TF

The Secure Tether provides high security mounting of Caregiver® with a 12-foot extended length coiled plastic-coated steel cable with security fittings. One end of the tether is securely attached to the Caregiver thermometer. The other end can be attached to the Wall Mount WMT-TF or to user's equipment. The included SecureMount $^{\text{TM}}$ key assures that attachment and removal is only by personnel holding the key.



12-FT. COILED CORD MODEL CC12-TF

Strong, flexible and plastic coated for easy disinfection/cleaning. The slip ring at the end of the coiled cord can be attached to the security pin on the bottom of the wall mount. The opposite end is attached to the fitting on the Caregiver thermometer.

